

Synergy^{RF}™ Bipolar Ablation System With Apollo^{RF}® MP90 Probe

Minimally Invasive and Effective Tissue Ablation and
Coagulation in Arthroscopic Procedures




Arthrex®

Synergy^{RF}™ Bipolar Ablation System With Apollo^{RF}® MP90 Probe

Features and Benefits:

- 90° tip configured to access and ablate tissue rapidly
- Easy-to-use, ergonomic hand controls designed to provide a seamless ablation experience
- Controlled, plasma-based bipolar RF technology
- Contains a durable one-piece active electrode
- Multiple ports for efficient bubble removal



Apollo^{RF} MP90 Probe



Synergy^{RF} Bipolar Ablation System



Synergy^{RF} Console (AR-9800):

- Integrates with Synergy Imaging System or operates as a stand-alone unit
- Probe-recognition function identifies probe currently in use
- Metal proximity detection function suspends power when in close proximity to scope and metal instruments
- User-friendly LCD touchscreen interface
- Easy-to-adjust ablation and coagulation settings



Synergy^{RF} Footswitch (AR-9800-F):

- Separate ablation and coagulation pedals
- Central ablation-level-adjustment button
- Integrated carrying handle

Ordering Information

Product Description	Item Number
Apollo ^{RF} MP90 probe	AR-9811
Synergy ^{RF} console	AR-9800
Synergy ^{RF} footswitch	AR-9800-F



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



Arthrex manufacturer,
authorized representative,
and importer information
(Arthrex eFUs)



US patent
information

arthrex.com