


Synergy OR Command™ System v2.0

OR Efficiency: Simplified




Synergy OR Command™ System v2.0

Arthrex is pleased to announce the latest release of the Synergy OR Command system v2.0. This release includes a new streaming engine designed to improve performance and reliability, and several new features to enhance usability.

A powerful tool for maximizing OR time, Synergy OR Command room status provides surgeons and staff with a centralized dashboard view of simultaneous

video streams across all operating rooms. Authorized users can access Synergy OR Command room status through a networked browser or smart device. Benefits of Synergy OR Command room status include adaptive scheduling, minimized downtime, accelerated surgeon workflows, and real-time information for decision-makers. Whether accessed on a computer, tablet, or wall display in an OR corridor, OR Command enhances efficiency in today's busy ORs.



Key Benefits of Synergy OR Command Room Status

- **New:** Updated streaming engine offers improved performance and reliability
- **New:** Manually arrange stream windows in Flat View
- **New:** Privacy Mode now available for each individual source
- **New:** Encrypted video streams provide enhanced cybersecurity
- **New:** 3x4 dashboard layout shows up to 12 simultaneous feeds
- Improve workflow efficiencies with OR visibility
- Kiosk Mode for dedicated displays
- Configurable dashboard layout options
- Scope view from Synergy camera consoles

Update Path for Existing OR Command System v1.0 Customers

Existing OR Command system v1.0 customers can update their system to the latest release in one of three ways:

1. Purchasing a Synergy OR Command support contract
2. Purchasing a Synergy Vision™ console and streaming license (included in the Workflow Support module and Premium software bundle)
3. Exercising their 1-year warranty, which is included in every Synergy OR Command room status installation

Ordering Information

Synergy OR Command room status software license	AR-3200-0044
Synergy OR Command room status support plan	SV-HWSW1-4OR Command1Y