

Variable Angle Locking (VAL) KreuLock™ Screws

Patient Information Leaflet



Helping Surgeons Treat Their Patients Better[®]

Since its inception, Arthrex has been committed to one mission: Helping Surgeons Treat Their Patients Better. We are strategically focused on constant product innovation through scientific research, surgeon collaboration, and medical education to make less invasive surgical procedures simpler, safer, and more reproducible. Each year, we develop more than 1000 new innovative products and procedures to advance minimally invasive orthopedics worldwide.

Arthrex has always remained a privately held company, which allows for the rapid evaluation of new technologies and ideas. Our economic strength enables us to develop products and techniques that truly make a difference without compromising on quality. Our experienced team of dedicated professionals represents a shared passion and commitment to delivering uncompromising quality to the health care providers who use our products and the millions of patients whose lives we impact.

The medical significance of our contributions serves as our primary benchmark of success and will continue into the future as the legacy of Arthrex.

Table of Contents

04	Device Description
05	Material Specifications
06	Postoperative Care
06	MRI Safety Information
07	Arthrex Standalone Screw Device Family
08	Screws with Plates used in Foot and Ankle Region
09	Screws used with the Arthrex Variable Angle (VA) Proximal Tibia Plating System
10	Screws with Plates used in the Arm and Elbow Region
11	Screws with Plates used in Hand, Wrist, and Humerus Region
12	Screws with Plates used in Clavicle Region
12	Implant Models
13	VAL KreuLock™ Compression Screws (2.0-3.0 mm solid)
15	VAL KreuLock™ Screw, 2.7 mm, Stainless Steel (SS) series
17	VAL KreuLock™ Screw, 2.7 mm, Titanium (Ti) series
19	VAL KreuLock Screw, 3.0 mm, Hybrid, Titanium (Ti) series
21	VAL KreuLock Screw, 3.0 mm, Titanium (Ti) series
23	VAL KreuLock™ Screw, 3.5 mm, Stainless Steel (SS) series
26	VAL KreuLock™ Screw, 3.5 mm, Titanium (Ti) series, Reinforced (Ri)
29	VAL KreuLock™ Screw, 3.5 mm, Titanium (Ti) series
32	Contact Information

Device Description

The Variable Angle Locking (VAL) KreuLock Screws consist of permanent implant devices manufactured from titanium or stainless steel alloy. The screw devices are available as solid or cannulated implants that are available in a variety of sizes. The screws are sold single-use and non-sterile or sterile.

The VAL KreuLock™ Compression Screws (2.0-3.0 mm solid)

The 2.0-3.0 mm solid VAL KreuLock™ Compression Screws are intended to be used in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, tibia, fibula, femur, pelvis, acetabulum, metacarpals, metatarsals, humerus, radius, ulna, calcaneus, clavicle and scapula. The screw may be used with the Arthrex Low Profile Plate, Small Fragment Plates, Distal Extremity Plates, Mesh Plates (2.7-3.0) and Distal Radius Plates.

The VAL KreuLock™ Compression Screws (2.7 mm solid)

The 2.7 mm solid VAL KreuLock Compression Screws are intended to be used in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, tibia, fibula, femur, pelvis, acetabulum, metacarpals, metatarsals, humerus, radius, ulna, calcaneus, and clavicle. When used with a plate, the screw may be used with the Arthrex Low Profile Plate, Small Fragment Plates, Distal Extremity Plates, Mesh Plates (2.7) and Distal Radius Plates.

The VAL KreuLock™ Compression Screws (3.5 mm and larger, solid)

The 3.5 mm and larger, solid VAL KreuLock Compression Screws are intended to be used in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, pelvis, acetabulum, metacarpals, metatarsals, femur, and fibula. The screws may be used with the Arthrex Low Profile Plates, Small Fragment Plates, Fracture Plates, Distal Extremity Plates, Distal Radius Plates, Humeral Fracture Plates, Osteotomy Plates, and Ankle Fusion Plates.

The Hybrid KreuLock™ Compression Screws (3.5 mm and larger, solid)

The 3.5 mm and larger, solid Hybrid KreuLock™ Compression Screws are intended to be used in a plate-screw system for internal bone fixation for bone fractures, fusions, Osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, pelvis, acetabulum, metacarpals, metatarsals, femur and fibula. The screws may be used with the Arthrex Low Profile Plate, Small Fragment Plates, Fracture Plates, Distal Extremity Plates, Distal Radius Plates, Humeral Fracture Plates, Osteotomy Plates, and Ankle Fusion Plates.

Material Specifications

The VAL KreuLock Screws are manufactured from titanium or stainless steel alloy.

Screws that are manufactured from titanium alloy (ASTM F136) contain:

- Titanium, (balance)
- Aluminum, (5.5 - 6.5%)
- Vanadium, (3.5 - 4.5%)
- Iron, (0.25% max)
- Oxygen, (0.13% max)
- *Other materials may be present at trace levels.

Screws that are manufactured from Stainless Steel, 316L (ASTM F 138) contain:

- Iron, (balance)
- Chromium, (17.0 - 19.0%)
- Nickel, (13.0 - 15.0%)
- Molybdenum, (2.25 - 3.00%)
- Manganese, (2.0%)
- Silicon, (1.0%)
- *Other materials may be present at trace levels.

Postoperative Care

Postoperative management is patient-specific and dependent on your doctor's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. Please be aware that surgery and recovery protocol may vary for each individual and any questions pertaining to the surgical procedure or postoperative protocol should be discussed with your surgeon.

Please call your doctor if:

- › You experience loss of function/range of motion
- › You develop a fever greater than 38 °C (100.4 °F)
- › Drainage continues from the site of your incision
- › Your surgical site becomes more swollen, tender, and painful, with increased difficulty performing your exercises.

If you have difficulty breathing or develop severe pain or chest pain, call 911 or report immediately to your local emergency room.

MRI Safety Information



Arthrex Standalone Screw Device Family

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when receiving an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with Arthrex Implants may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Family	Arthrex Standalone Screw Device Family
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3000 gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR (Specific Absorption Rate)	2 W/kg
Maximum Head SAR	3.2 W/kg
Scan Regions	The edge of the device should be positioned at least 35 cm away from the isocenter
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF.
MR Image Artifact	The presence of this implant may produce an image artifact.

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

If information about a specific parameter is not included, there are no conditions associated with that parameter.



A person with Arthrex Standalone Screws can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use (<https://edfu.arthrex.com>) or by calling Arthrex customer service at ☎ +1 800 934-4404.

MRI Safety Information



Screws with Plates used in Foot and Ankle Region

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when receiving an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with Arthrex Implants may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Name	Screws with Plates used in Foot and Ankle region
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3,000 gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Maximum Whole-Body (WB) Specific Absorption Rate (SAR) when imaging landmark is < 35 cm from the implant device	1 W/kg with a scan time restriction of 10 minutes, followed by rest of 20 minutes for a total scan time of 60 minutes (three 10-minute scans)*
Maximum WB SAR when imaging landmark is ≥35 cm from the implant device	2 W/kg with 60 minutes of continuous scanning
Maximum Head SAR	3.2 W/kg**
MR Image Artifact	The presence of this implant may produce an image artifact.

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

If information about a specific parameter is not included, there are no conditions associated with that parameter.

*Actual scanning time (continuous RF exposure) is still less than 60-minutes

**Device is not located near the head region. Therefore, head SAR region is not applicable and can be labeled at 3.2 W/kg without any safety concern.



A person with Screws with Plates used in Foot and Ankle region can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use

(<https://edfu.arthrex.com>) or by calling Arthrex customer service at ☎ +1 800 934-4404.

MRI Safety Information



Screws used with the Arthrex Variable Angle (VA) Proximal Tibia Plating System

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when getting an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with Arthrex Implants may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Name	Screws used with the Arthrex Variable Angle (VA) Proximal Tibia Plating System
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3000 gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body (WB) Specific Absorption Rate (SAR) when imaging landmark is < 35 cm from the implant device	For 1.5T: 2 W/kg with 60 minutes of continuous scanning For 3T: 1 W/kg with 30 minutes of continuous scanning, followed by a rest period of 15 minutes. This scan-rest sequence may be repeated once.
Maximum WB SAR when imaging landmark is ≥35 cm from the implant device	2 W/kg with 60 minutes of continuous scanning
Maximum Head SAR	3.2 W/kg*
MR Image Artifact	The presence of this implant may produce an image artifact.
Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.	
If information about a specific parameter is not included, there are no conditions associated with that parameter.	
*Device is not located near the head region. Therefore, head SAR region is not applicable and can be labeled at 3.2 W/kg without any safety concern.	



A person with Screws used with the Arthrex VA Proximal Tibia Plating System can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use (<https://edfu.arthrex.com>) or by calling Arthrex customer service at ☎ +1 800 934-4404.

MRI Safety Information

Screws with Plates used in the Arm and Elbow Region

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when getting an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with Arthrex Implants may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Name	Screws with Plates used in Arm and Elbow region
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3000 Gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Maximum Whole-Body SAR (Specific Absorption Rate)	2 W/kg
Maximum Head SAR	N/A
Scan Duration	Under the defined scan conditions, the Arthrex elbow plates can be scanned continuously for 15 minutes. Following a 15-minute cooling time, the sequence can be repeated twice within 60 minutes.
MR Image Artifact	The presence of this implant may produce an image artifact.

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

If information about a specific parameter is not included, there are no conditions associated with that parameter.



A person with Screws with Plates used in the Arm and Elbow region can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use (<https://edfu.arthrex.com>) or by calling Arthrex customer service at +1 800 934-4404.

MRI Safety Information

Screws with Plates used in Hand, Wrist, and Humerus Region

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when getting an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with Arthrex Implants may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Family	Screws with Plates used in Hand, Wrist, and Humerus region
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3000 Gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR (Specific Absorption Rate)	1 W/kg (Normal Operating Mode)
Maximum Head SAR	N/A

Scan Duration	Under the scan conditions defined above, the patients can be scanned continuously for 15 minutes with a 5-minute wait time. The sequence can be repeated three times within 60 minutes.
---------------	---

MR Image Artifact	The presence of this implant may produce an image artifact.
-------------------	---

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

If information about a specific parameter is not included, there are no conditions associated with that parameter.



A person with Screws with Plates used in Hand, Wrist, and Humerus region can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use (<https://edfu.arthrex.com>) or by calling Arthrex customer service at ☎ +1 800 934-4404.

MRI Safety Information

Screws with Plates used in Clavicle Region

MRI, or Magnetic Resonance Imaging, is an imaging technique utilizing a strong magnetic field to produce detailed anatomical images. This section details the information that you should be aware of when getting an MRI scan.

1. MR (Magnetic Resonance) Conditional



MRI Safety Information

A person with the Arthrex Clavicle Plating system may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Family	Screws with Plates used in Clavicle region
Static Magnetic Field Strength (B ₀)	1.5 Tesla (T) or 3.0 T
Maximum Spatial Field Gradient	30 T/m or 3000 Gauss/cm
RF (Radio Frequency) Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Volume RF body coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR (Specific Absorption Rate)	2 W/kg (normal operating mode)
Maximum Head SAR	N/A
Scan Duration	Under the scan conditions defined, the Arthrex Clavicle Plates can be scanned continuously for 60 minutes.
MR Image Artifact	The presence of this implant may produce an image artifact.

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. Do not conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

If information about a specific parameter is not included, there are no conditions associated with that parameter.



A person with Screws with Plates used in Clavicle region can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in severe injury. Full MRI safety information is available in the MRI Safety Information section of this patient information leaflet, Directions for Use

(<https://edfu.arthrex.com>) or by calling Arthrex customer service at ☎ +1 800 934-4404.

Implant Models

VAL KreuLock™ Compression Screws (2.0-3.0 mm solid)

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Part Number	Part Description
AR-18720VCL-06	VAL KreuLock™ Screw, Ti, 2.0 mm x 6 mm
AR-18720VCL-06S	VAL KreuLock™ Screw, Ti, 2.0 mm x 6 mm, Strl
AR-18720VCL-07	VAL KreuLock™ Screw, Ti, 2.0 mm x 7 mm
AR-18720VCL-07S	VAL KreuLock™ Screw, Ti, 2.0 mm x 7 mm, Strl
AR-18720VCL-08	VAL KreuLock™ Screw, Ti, 2.0 mm x 8 mm
AR-18720VCL-08S	VAL KreuLock™ Screw, Ti, 2.0 mm x 8 mm, Strl
AR-18720VCL-09	VAL KreuLock™ Screw, Ti, 2.0 mm x 9 mm
AR-18720VCL-09S	VAL KreuLock™ Screw, Ti, 2.0 mm x 9 mm, Strl
AR-18720VCL-10	VAL KreuLock™ Screw, Ti, 2.0 mm x 10 mm
AR-18720VCL-10S	VAL KreuLock™ Screw, Ti, 2.0 mm x 10 mm, Strl
AR-18720VCL-11	VAL KreuLock™ Screw, Ti, 2.0 mm x 11 mm
AR-18720VCL-11S	VAL KreuLock™ Screw, Ti, 2.0 mm x 11 mm, Strl
AR-18720VCL-12	VAL KreuLock™ Screw, Ti, 2.0 mm x 12 mm
AR-18720VCL-12S	VAL KreuLock™ Screw, Ti, 2.0 mm x 12 mm, Strl
AR-18720VCL-13	VAL KreuLock™ Screw, Ti, 2.0 mm x 13 mm
AR-18720VCL-13S	VAL KreuLock™ Screw, Ti, 2.0 mm x 13 mm, Strl
AR-18720VCL-14	VAL KreuLock™ Screw, Ti, 2.0 mm x 14 mm
AR-18720VCL-14S	VAL KreuLock™ Screw, Ti, 2.0 mm x 14 mm, Strl
AR-18720VCL-15	VAL KreuLock™ Screw, Ti, 2.0 mm x 15 mm
AR-18720VCL-15S	VAL KreuLock™ Screw, Ti, 2.0 mm x 15 mm, Strl
AR-18720VCL-16	VAL KreuLock™ Screw, Ti, 2.0 mm x 16 mm
AR-18720VCL-16S	VAL KreuLock™ Screw, Ti, 2.0 mm x 16 mm, Strl
AR-18720VCL-18	VAL KreuLock™ Screw, Ti, 2.0 mm x 18 mm
AR-18720VCL-18S	VAL KreuLock™ Screw, Ti, 2.0 mm x 18 mm, Strl
AR-18720VCL-20	VAL KreuLock™ Screw, Ti, 2.0 mm x 20 mm
AR-18720VCL-20S	VAL KreuLock™ Screw, Ti, 2.0 mm x 20 mm, Strl
AR-18720VCL-22	VAL KreuLock™ Screw, Ti, 2.0 mm x 22 mm
AR-18720VCL-22S	VAL KreuLock™ Screw, Ti, 2.0 mm x 22 mm, Strl
AR-18720VCL-24	VAL KreuLock™ Screw, Ti, 2.0 mm x 24 mm
AR-18720VCL-24S	VAL KreuLock™ Screw, Ti, 2.0 mm x 24 mm, Strl
AR-18720VCL-26	VAL KreuLock™ Screw, Ti, 2.0 mm x 26 mm
AR-18720VCL-26S	VAL KreuLock™ Screw, Ti, 2.0 mm x 26 mm, Strl
AR-18720VCL-28	VAL KreuLock™ Screw, Ti, 2.0 mm x 28 mm
AR-18720VCL-28S	VAL KreuLock™ Screw, Ti, 2.0 mm x 28 mm, Strl

Part Number	Part Description
AR-18720VCL-30	VAL KreuLock™ Screw, Ti, 2.0 mm x 30 mm
AR-18720VCL-30S	VAL KreuLock™ Screw, Ti, 2.0 mm x 30 mm, Strl
AR-18720VCL-32	VAL KreuLock™ Screw, Ti, 2.0 mm x 32 mm
AR-18720VCL-32S	VAL KreuLock™ Screw, Ti, 2.0 mm x 32 mm, Strl
AR-18720VCL-34	VAL KreuLock™ Screw, Ti, 2.0 mm x 34 mm
AR-18720VCL-34S	VAL KreuLock™ Screw, Ti, 2.0 mm x 34 mm, Strl
AR-18720VCL-36	VAL KreuLock™ Screw, Ti, 2.0 mm x 36 mm
AR-18720VCL-36S	VAL KreuLock™ Screw, Ti, 2.0 mm x 36 mm, Strl
AR-18720VCL-38	VAL KreuLock™ Screw, Ti, 2.0 mm x 38 mm
AR-18720VCL-38S	VAL KreuLock™ Screw, Ti, 2.0 mm x 38 mm, Strl
AR-18720VCL-40	VAL KreuLock™ Screw, Ti, 2.0 mm x 40 mm
AR-18720VCL-40S	VAL KreuLock™ Screw, Ti, 2.0 mm x 40 mm, Strl

Implant Models

VAL KreuLock™ Screw, 2.7 mm, Stainless Steel (SS) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8827VCL-10	VAL KreuLock™ Comp. Screw, SS, 2.7 x 10 mm
AR-8827VCL-10S	VAL KreuLock™ Screw, SS, 2.7 x 10 mm, Strl
AR-8827VCL-11	VAL KreuLock™ Comp. Screw, SS, 2.7 x 11 mm
AR-8827VCL-11S	VAL KreuLock™ Screw, SS, 2.7 x 11 mm, Strl
AR-8827VCL-12	VAL KreuLock™ Comp. Screw, SS, 2.7 x 12 mm
AR-8827VCL-12S	VAL KreuLock™ Screw, SS, 2.7 x 12 mm, Strl
AR-8827VCL-13	VAL KreuLock™ Comp. Screw, SS, 2.7 x 13 mm
AR-8827VCL-13S	VAL KreuLock™ Screw, SS, 2.7 x 13 mm, Strl
AR-8827VCL-14	VAL KreuLock™ Comp. Screw, SS, 2.7 x 14 mm
AR-8827VCL-14S	VAL KreuLock™ Screw, SS, 2.7 x 14 mm, Strl
AR-8827VCL-15	VAL KreuLock™ Comp. Screw, SS, 2.7 x 15 mm
AR-8827VCL-15S	VAL KreuLock™ Screw, SS, 2.7 x 15 mm, Strl
AR-8827VCL-16	VAL KreuLock™ Comp. Screw, SS, 2.7 x 16 mm
AR-8827VCL-16S	VAL KreuLock™ Screw, SS, 2.7 x 16 mm, Strl
AR-8827VCL-18	VAL KreuLock™ Comp. Screw, SS, 2.7 x 18 mm
AR-8827VCL-18S	VAL KreuLock™ Screw, SS, 2.7 x 18 mm, Strl
AR-8827VCL-20	VAL KreuLock™ Comp. Screw, SS, 2.7 x 20 mm
AR-8827VCL-20S	VAL KreuLock™ Screw, SS, 2.7 x 20 mm, Strl
AR-8827VCL-22	VAL KreuLock™ Comp. Screw, SS, 2.7 x 22 mm
AR-8827VCL-22S	VAL KreuLock™ Screw, SS, 2.7 x 22 mm, Strl
AR-8827VCL-24	VAL KreuLock™ Comp. Screw, SS, 2.7 x 24 mm
AR-8827VCL-24S	VAL KreuLock™ Screw, SS, 2.7 x 24 mm, Strl
AR-8827VCL-26	VAL KreuLock™ Comp. Screw, SS, 2.7 x 26 mm
AR-8827VCL-26S	VAL KreuLock™ Screw, SS, 2.7 x 26 mm, Strl
AR-8827VCL-28	VAL KreuLock™ Comp. Screw, SS, 2.7 x 28 mm
AR-8827VCL-28S	VAL KreuLock™ Screw, SS, 2.7 x 28 mm, Strl
AR-8827VCL-30	VAL KreuLock™ Comp. Screw, SS, 2.7 x 30 mm
AR-8827VCL-30S	VAL KreuLock™ Screw, SS, 2.7 x 30 mm, Strl
AR-8827VCL-32	VAL KreuLock™ Comp. Screw, SS, 2.7 x 32 mm
AR-8827VCL-32S	VAL KreuLock™ Screw, SS, 2.7 x 32 mm, Strl
AR-8827VCL-34	VAL KreuLock™ Comp. Screw, SS, 2.7 x 34 mm
AR-8827VCL-34S	VAL KreuLock™ Screw, SS, 2.7 x 34 mm, Strl
AR-8827VCL-36	VAL KreuLock™ Comp. Screw, SS, 2.7 x 36 mm
AR-8827VCL-36S	VAL KreuLock™ Screw, SS, 2.7 x 36 mm, Strl

Item Number	Product Description
AR-8827VCL-38	VAL KreuLock™ Comp. Screw, SS, 2.7 x 38 mm
AR-8827VCL-38S	VAL KreuLock™ Screw, SS, 2.7 x 38 mm, Strl
AR-8827VCL-40	VAL KreuLock™ Comp. Screw, SS, 2.7 x 40 mm
AR-8827VCL-40S	VAL KreuLock™ Screw, SS, 2.7 x 40 mm, Strl
AR-8827VCL-42	VAL KreuLock™ Comp. Screw, SS, 2.7 x 42 mm
AR-8827VCL-42S	VAL KreuLock™ Screw, SS, 2.7 x 42 mm, Strl
AR-8827VCL-44	VAL KreuLock™ Comp. Screw, SS, 2.7 x 44 mm
AR-8827VCL-44S	VAL KreuLock™ Screw, SS, 2.7 x 44 mm, Strl
AR-8827VCL-45	VAL KreuLock™ Comp. Screw, SS, 2.7 x 45 mm
AR-8827VCL-45S	VAL KreuLock™ Screw, SS, 2.7 x 45 mm, Strl
AR-8827VCL-46	VAL KreuLock™ Comp. Screw, SS, 2.7 x 46 mm
AR-8827VCL-46S	VAL KreuLock™ Screw, SS, 2.7 x 46 mm, Strl
AR-8827VCL-48	VAL KreuLock™ Comp. Screw, SS, 2.7 x 48 mm
AR-8827VCL-48S	VAL KreuLock™ Screw, SS, 2.7 x 48 mm, Strl
AR-8827VCL-50	VAL KreuLock™ Comp. Screw, SS, 2.7 x 50 mm
AR-8827VCL-50S	VAL KreuLock™ Screw, SS, 2.7 x 50 mm, Strl
AR-8827VCL-52	VAL KreuLock™ Comp. Screw, SS, 2.7 x 52 mm
AR-8827VCL-52S	VAL KreuLock™ Screw, SS, 2.7 x 52 mm, Strl
AR-8827VCL-54	VAL KreuLock™ Comp. Screw, SS, 2.7 x 54 mm
AR-8827VCL-54S	VAL KreuLock™ Screw, SS, 2.7 x 54 mm, Strl
AR-8827VCL-55	VAL KreuLock™ Comp. Screw, SS, 2.7 x 55 mm
AR-8827VCL-55S	VAL KreuLock™ Screw, SS, 2.7 x 55 mm, Strl
AR-8827VCL-56	VAL KreuLock™ Comp. Screw, SS, 2.7 x 56 mm
AR-8827VCL-56S	VAL KreuLock™ Screw, SS, 2.7 x 56 mm, Strl
AR-8827VCL-58	VAL KreuLock™ Comp. Screw, SS, 2.7 x 58 mm
AR-8827VCL-58S	VAL KreuLock™ Screw, SS, 2.7 x 58 mm, Strl
AR-8827VCL-60	VAL KreuLock™ Comp. Screw, SS, 2.7 x 60 mm
AR-8827VCL-60S	VAL KreuLock™ Screw, SS, 2.7 x 60 mm, Strl
AR-8827VCL-65	VAL KreuLock™ Comp. Screw, SS, 2.7 x 65 mm
AR-8827VCL-65S	VAL KreuLock™ Screw, SS, 2.7 x 65 mm, Strl
AR-8827VCL-70	VAL KreuLock™ Comp. Screw, SS, 2.7 x 70 mm
AR-8827VCL-70S	VAL KreuLock™ Screw, SS, 2.7 x 70 mm, Strl
AR-8827VCL-75	VAL KreuLock™ Comp. Screw, SS, 2.7 x 75 mm
AR-8827VCL-75S	VAL KreuLock™ Screw, SS, 2.7 x 75 mm, Strl
AR-8827VCL-80	VAL KreuLock™ Comp. Screw, SS, 2.7 x 80 mm
AR-8827VCL-80S	VAL KreuLock™ Screw, SS, 2.7 x 80 mm, Strl

Implant Models

VAL KreuLock™ Screw, 2.7 mm, Titanium (Ti) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8927VCL-10	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 10 mm
AR-8927VCL-10S	VAL KreuLock™ Screw, Ti, 2.7 x 10 mm, Strl
AR-8927VCL-11	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 11 mm
AR-8927VCL-11S	VAL KreuLock™ Screw, Ti, 2.7 x 11 mm, Strl
AR-8927VCL-12	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 12 mm
AR-8927VCL-12S	VAL KreuLock™ Screw, Ti, 2.7 x 12 mm, Strl
AR-8927VCL-13	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 13 mm
AR-8927VCL-13S	VAL KreuLock™ Screw, Ti, 2.7 x 13 mm, Strl
AR-8927VCL-14	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 14 mm
AR-8927VCL-14S	VAL KreuLock™ Screw, Ti, 2.7 x 14 mm, Strl
AR-8927VCL-15	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 15 mm
AR-8927VCL-15S	VAL KreuLock™ Screw, Ti, 2.7 x 15 mm, Strl
AR-8927VCL-16	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 16 mm
AR-8927VCL-16S	VAL KreuLock™ Screw, Ti, 2.7 x 16 mm, Strl
AR-8927VCL-18	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 18 mm
AR-8927VCL-18S	VAL KreuLock™ Screw, Ti, 2.7 x 18 mm, Strl
AR-8927VCL-20	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 20 mm
AR-8927VCL-20S	VAL KreuLock™ Screw, Ti, 2.7 x 20 mm, Strl
AR-8927VCL-22	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 22 mm
AR-8927VCL-22S	VAL KreuLock™ Screw, Ti, 2.7 x 22 mm, Strl
AR-8927VCL-24	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 24 mm
AR-8927VCL-24S	VAL KreuLock™ Screw, Ti, 2.7 x 24 mm, Strl
AR-8927VCL-26	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 26 mm
AR-8927VCL-26S	VAL KreuLock™ Screw, Ti, 2.7 x 26 mm, Strl
AR-8927VCL-28	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 28 mm
AR-8927VCL-28S	VAL KreuLock™ Screw, Ti, 2.7 x 28 mm, Strl
AR-8927VCL-30	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 30 mm
AR-8927VCL-30S	VAL KreuLock™ Screw, Ti, 2.7 x 30 mm, Strl
AR-8927VCL-32	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 32 mm
AR-8927VCL-32S	VAL KreuLock™ Screw, Ti, 2.7 x 32 mm, Strl
AR-8927VCL-34	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 34 mm
AR-8927VCL-34S	VAL KreuLock™ Screw, Ti, 2.7 x 34 mm, Strl
AR-8927VCL-36	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 36 mm
AR-8927VCL-36S	VAL KreuLock™ Screw, Ti, 2.7 x 36 mm, Strl

Item Number	Product Description
AR-8927VCL-38	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 38 mm
AR-8927VCL-38S	VAL KreuLock™ Screw, Ti, 2.7 x 38 mm, Strl
AR-8927VCL-40	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 40 mm
AR-8927VCL-40S	VAL KreuLock™ Screw, Ti, 2.7 x 40 mm, Strl
AR-8927VCL-42	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 42 mm
AR-8927VCL-42S	VAL KreuLock™ Screw, Ti, 2.7 x 42 mm, Strl
AR-8927VCL-44	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 44 mm
AR-8927VCL-44S	VAL KreuLock™ Screw, Ti, 2.7 x 44 mm, Strl
AR-8927VCL-45	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 45 mm
AR-8927VCL-45S	VAL KreuLock™ Screw, Ti, 2.7 x 45 mm, Strl
AR-8927VCL-46	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 46 mm
AR-8927VCL-46S	VAL KreuLock™ Screw, Ti, 2.7 x 46 mm, Strl
AR-8927VCL-48	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 48 mm
AR-8927VCL-48S	VAL KreuLock™ Screw, Ti, 2.7 x 48 mm, Strl
AR-8927VCL-50	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 50 mm
AR-8927VCL-50S	VAL KreuLock™ Screw, Ti, 2.7 x 50 mm, Strl
AR-8927VCL-52	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 52 mm
AR-8927VCL-52S	VAL KreuLock™ Screw, Ti, 2.7 x 52 mm, Strl
AR-8927VCL-54	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 54 mm
AR-8927VCL-54S	VAL KreuLock™ Screw, Ti, 2.7 x 54 mm, Strl
AR-8927VCL-55	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 55 mm
AR-8927VCL-55S	VAL KreuLock™ Screw, Ti, 2.7 x 55 mm, Strl
AR-8927VCL-56	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 56 mm
AR-8927VCL-56S	VAL KreuLock™ Screw, Ti, 2.7 x 56 mm, Strl
AR-8927VCL-58	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 58 mm
AR-8927VCL-58S	VAL KreuLock™ Screw, Ti, 2.7 x 58 mm, Strl
AR-8927VCL-60	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 60 mm
AR-8927VCL-60S	VAL KreuLock™ Screw, Ti, 2.7 x 60 mm, Strl
AR-8927VCL-65	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 65 mm
AR-8927VCL-65S	VAL KreuLock™ Screw, Ti, 2.7 x 65 mm, Strl
AR-8927VCL-70	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 70 mm
AR-8927VCL-70S	VAL KreuLock™ Screw, Ti, 2.7 x 70 mm, Strl
AR-8927VCL-75	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 75 mm
AR-8927VCL-75S	VAL KreuLock™ Screw, Ti, 2.7 x 75 mm, Strl
AR-8927VCL-80	VAL KreuLock™ Comp. Screw, Ti, 2.7 x 80 mm
AR-8927VCL-80S	VAL KreuLock™ Screw, Ti, 2.7 x 80 mm, Strl

Implant Models

VAL KreuLock Screw, 3.0 mm, Hybrid, Titanium (Ti) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8933HVCL-12	Hybrid VAL KreuLock Screw, Ti, 3.0 x 12 mm
AR-8933HVCL-12S	Hyb VAL KreuLock Screw, Ti, 3.0 x 12 mm, Strl
AR-8933HVCL-14	Hybrid VAL KreuLock Screw, Ti, 3.0 x 14 mm
AR-8933HVCL-14S	Hyb VAL KreuLock Screw, Ti, 3.0 x 14 mm, Strl
AR-8933HVCL-16	Hybrid VAL KreuLock Screw, Ti, 3.0 x 16 mm
AR-8933HVCL-16S	Hyb VAL KreuLock Screw, Ti, 3.0 x 16 mm, Strl
AR-8933HVCL-18	Hybrid VAL KreuLock Screw, Ti, 3.0 x 18 mm
AR-8933HVCL-18S	Hyb VAL KreuLock Screw, Ti, 3.0 x 18 mm, Strl
AR-8933HVCL-20	Hybrid VAL KreuLock Screw, Ti, 3.0 x 20 mm
AR-8933HVCL-20S	Hyb VAL KreuLock Screw, Ti, 3.0 x 20 mm, Strl
AR-8933HVCL-22	Hybrid VAL KreuLock Screw, Ti, 3.0 x 22 mm
AR-8933HVCL-22S	Hyb VAL KreuLock Screw, Ti, 3.0 x 22 mm, Strl
AR-8933HVCL-24	Hybrid VAL KreuLock Screw, Ti, 3.0 x 24 mm
AR-8933HVCL-24S	Hyb VAL KreuLock Screw, Ti, 3.0 x 24 mm, Strl
AR-8933HVCL-26	Hybrid VAL KreuLock Screw, Ti, 3.0 x 26 mm
AR-8933HVCL-26S	Hyb VAL KreuLock Screw, Ti, 3.0 x 26 mm, Strl
AR-8933HVCL-28	Hybrid VAL KreuLock Screw, Ti, 3.0 x 28 mm
AR-8933HVCL-28S	Hyb VAL KreuLock Screw, Ti, 3.0 x 28 mm, Strl
AR-8933HVCL-30	Hybrid VAL KreuLock Screw, Ti, 3.0 x 30 mm
AR-8933HVCL-30S	Hyb VAL KreuLock Screw, Ti, 3.0 x 30 mm, Strl
AR-8933HVCL-32	Hybrid VAL KreuLock Screw, Ti, 3.0 x 32 mm
AR-8933HVCL-32S	Hyb VAL KreuLock Screw, Ti, 3.0 x 32 mm, Strl
AR-8933HVCL-34	Hybrid VAL KreuLock Screw, Ti, 3.0 x 34 mm
AR-8933HVCL-34S	Hyb VAL KreuLock Screw, Ti, 3.0 x 34 mm, Strl
AR-8933HVCL-36	Hybrid VAL KreuLock Screw, Ti, 3.0 x 36 mm
AR-8933HVCL-36S	Hyb VAL KreuLock Screw, Ti, 3.0 x 36 mm, Strl
AR-8933HVCL-38	Hybrid VAL KreuLock Screw, Ti, 3.0 x 38 mm
AR-8933HVCL-38S	Hyb VAL KreuLock Screw, Ti, 3.0 x 38 mm, Strl
AR-8933HVCL-40	Hybrid VAL KreuLock Screw, Ti, 3.0 x 40 mm
AR-8933HVCL-40S	Hyb VAL KreuLock Screw, Ti, 3.0 x 40 mm, Strl
AR-8933HVCL-42	Hybrid VAL KreuLock Screw, Ti, 3.0 x 42 mm
AR-8933HVCL-42S	Hyb VAL KreuLock Screw, Ti, 3.0 x 42 mm, Strl
AR-8933HVCL-44	Hybrid VAL KreuLock Screw, Ti, 3.0 x 44 mm
AR-8933HVCL-44S	Hyb VAL KreuLock Screw, Ti, 3.0 x 44 mm, Strl

Item Number	Product Description
AR-8933HVCL-46	Hybrid VAL KreuLock Screw, Ti, 3.0 x 46 mm
AR-8933HVCL-46S	Hyb VAL KreuLock Screw, Ti, 3.0 x 46 mm, Strl
AR-8933HVCL-48	Hybrid VAL KreuLock Screw, Ti, 3.0 x 48 mm
AR-8933HVCL-48S	Hyb VAL KreuLock Screw, Ti, 3.0 x 48 mm, Strl
AR-8933HVCL-50	Hybrid VAL KreuLock Screw, Ti, 3.0 x 50 mm
AR-8933HVCL-50S	Hyb VAL KreuLock Screw, Ti, 3.0 x 50 mm, Strl
AR-8933HVCL-52	Hybrid VAL KreuLock Screw, Ti, 3.0 x 52 mm
AR-8933HVCL-52S	Hyb VAL KreuLock Screw, Ti, 3.0 x 52 mm, Strl
AR-8933HVCL-54	Hybrid VAL KreuLock Screw, Ti, 3.0 x 54 mm
AR-8933HVCL-54S	Hyb VAL KreuLock Screw, Ti, 3.0 x 54 mm, Strl
AR-8933HVCL-56	Hybrid VAL KreuLock Screw, Ti, 3.0 x 56 mm
AR-8933HVCL-56S	Hyb VAL KreuLock Screw, Ti, 3.0 x 56 mm, Strl
AR-8933HVCL-58	Hybrid VAL KreuLock Screw, Ti, 3.0 x 58 mm
AR-8933HVCL-58S	Hyb VAL KreuLock Screw, Ti, 3.0 x 58 mm, Strl
AR-8933HVCL-60	Hybrid VAL KreuLock Screw, Ti, 3.0 x 60 mm
AR-8933HVCL-60S	Hyb VAL KreuLock Screw, Ti, 3.0 x 60 mm, Strl
AR-8933HVCL-65	Hyb VAL KreuLock Screw, 3.0 x 65 mm
AR-8933HVCL-65S	Hyb VAL KreuLock™ Screw, 3.0 x 65 mm, Strl
AR-8933HVCL-70	Hyb VAL KreuLock Screw, 3.0 x 70 mm
AR-8933HVCL-70S	Hyb VAL KreuLock™ Screw, 3.0 x 70 mm, Strl
AR-8933HVCL-75	Hyb VAL KreuLock Screw, 3.0 x 75 mm
AR-8933HVCL-75S	Hyb VAL KreuLock™ Screw, 3.0 x 75 mm, Strl
AR-8933HVCL-80	Hyb VAL KreuLock Screw, 3.0 x 80 mm
AR-8933HVCL-80S	Hyb VAL KreuLock™ Screw, 3.0 x 80 mm, Strl

Implant Models

VAL KreuLock Screw, 3.0 mm, Titanium (Ti) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8933VCL-10	VAL KreuLock Screw, Ti, 3.0 x 10 mm
AR-8933VCL-10S	VAL KreuLock Screw, Ti, 3.0 x 10 mm, Strl
AR-8933VCL-12	VAL KreuLock Screw, Ti, 3.0 x 12 mm
AR-8933VCL-12S	VAL KreuLock Screw, Ti, 3.0 x 12 mm, Strl
AR-8933VCL-14	VAL KreuLock Screw, Ti, 3.0 x 14 mm
AR-8933VCL-14S	VAL KreuLock Screw, Ti, 3.0 x 14 mm, Strl
AR-8933VCL-16	VAL KreuLock Screw, Ti, 3.0 x 16 mm
AR-8933VCL-16S	VAL KreuLock Screw, Ti, 3.0 x 16 mm, Strl
AR-8933VCL-18	VAL KreuLock Screw, Ti, 3.0 x 18 mm
AR-8933VCL-18S	VAL KreuLock Screw, Ti, 3.0 x 18 mm, Strl
AR-8933VCL-20	VAL KreuLock Screw, Ti, 3.0 x 20 mm
AR-8933VCL-20S	VAL KreuLock Screw, Ti, 3.0 x 20 mm, Strl
AR-8933VCL-22	VAL KreuLock Screw, Ti, 3.0 x 22 mm
AR-8933VCL-22S	VAL KreuLock Screw, Ti, 3.0 x 22 mm, Strl
AR-8933VCL-24	VAL KreuLock Screw, Ti, 3.0 x 24 mm
AR-8933VCL-24S	VAL KreuLock Screw, Ti, 3.0 x 24 mm, Strl
AR-8933VCL-26	VAL KreuLock Screw, Ti, 3.0 x 26 mm
AR-8933VCL-26S	VAL KreuLock Screw, Ti, 3.0 x 26 mm, Strl
AR-8933VCL-28	VAL KreuLock Screw, Ti, 3.0 x 28 mm
AR-8933VCL-28S	VAL KreuLock Screw, Ti, 3.0 x 28 mm, Strl
AR-8933VCL-30	VAL KreuLock Screw, Ti, 3.0 x 30 mm
AR-8933VCL-30S	VAL KreuLock Screw, Ti, 3.0 x 30 mm, Strl
AR-8933VCL-32	VAL KreuLock Screw, Ti, 3.0 x 32 mm
AR-8933VCL-32S	VAL KreuLock Screw, Ti, 3.0 x 32 mm, Strl
AR-8933VCL-34	VAL KreuLock Screw, Ti, 3.0 x 34 mm
AR-8933VCL-34S	VAL KreuLock Screw, Ti, 3.0 x 34 mm, Strl
AR-8933VCL-36	VAL KreuLock Screw, Ti, 3.0 x 36
AR-8933VCL-36S	VAL KreuLock Screw, Ti, 3.0 x 36, Strl
AR-8933VCL-38	VAL KreuLock Screw, Ti, 3.0 x 38 mm
AR-8933VCL-38S	VAL KreuLock Screw, Ti, 3.0 x 38 mm, Strl
AR-8933VCL-40	VAL KreuLock Screw, Ti, 3.0 x 40 mm
AR-8933VCL-40S	VAL KreuLock Screw, Ti, 3.0 x 40 mm, Strl
AR-8933VCL-42	VAL KreuLock Screw, Ti, 3.0 x 42 mm
AR-8933VCL-42S	VAL KreuLock™ Screw, Ti, 3.0 x 42 mm, Strl

Item Number	Product Description
AR-8933VCL-44	VAL KreuLock Screw, Ti, 3.0 x 44 mm
AR-8933VCL-44S	VAL KreuLock™ Screw, Ti, 3.0 x 44 mm, Strl
AR-8933VCL-46	VAL KreuLock Screw, Ti, 3.0 x 46 mm
AR-8933VCL-46S	VAL KreuLock™ Screw, Ti, 3.0 x 46 mm, Strl
AR-8933VCL-48	VAL KreuLock Screw, Ti, 3.0 x 48 mm
AR-8933VCL-48S	VAL KreuLock™ Screw, Ti, 3.0 x 48 mm, Strl
AR-8933VCL-50	VAL KreuLock Screw, Ti, 3.0 x 50 mm
AR-8933VCL-50S	VAL KreuLock™ Screw, Ti, 3.0 x 50 mm, Strl
AR-8933VCL-52	VAL KreuLock Screw, Ti, 3.0 x 52 mm
AR-8933VCL-52S	VAL KreuLock™ Screw, Ti, 3.0 x 52 mm, Strl
AR-8933VCL-54	VAL KreuLock Screw, Ti, 3.0 x 54 mm
AR-8933VCL-54S	VAL KreuLock™ Screw, Ti, 3.0 x 54 mm, Strl
AR-8933VCL-56	VAL KreuLock Screw, Ti, 3.0 x 56 mm
AR-8933VCL-56S	VAL KreuLock™ Screw, Ti, 3.0 x 56 mm, Strl
AR-8933VCL-58	VAL KreuLock Screw, Ti, 3.0 x 58 mm
AR-8933VCL-58S	VAL KreuLock™ Screw, Ti, 3.0 x 58 mm, Strl
AR-8933VCL-60	VAL KreuLock Screw, Ti, 3.0 x 60 mm
AR-8933VCL-60S	VAL KreuLock™ Screw, Ti, 3.0 x 60 mm, Strl
AR-8933VCL-65	VAL KreuLock Screw, Ti, 3.0 x 65 mm
AR-8933VCL-65S	VAL KreuLock™ Screw, Ti, 3.0 x 65 mm, Strl
AR-8933VCL-70	VAL KreuLock Screw, Ti, 3.0 x 70 mm
AR-8933VCL-70S	VAL KreuLock™ Screw, Ti, 3.0 x 70 mm, Strl
AR-8933VCL-75	VAL KreuLock Screw, Ti, 3.0 x 75 mm
AR-8933VCL-75S	VAL KreuLock™ Screw, Ti, 3.0 x 75 mm, Strl
AR-8933VCL-80	VAL KreuLock Screw, Ti, 3.0 x 80 mm
AR-8933VCL-80S	VAL KreuLock™ Screw, Ti, 3.0 x 80 mm, Strl

Implant Models

VAL KreuLock™ Screw, 3.5 mm, Stainless Steel (SS) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8835VCL-10	VAL KreuLock™ Comp. Screw, SS, 3.5 x 10 mm
AR-8835VCL-10S	VAL KreuLock™ Screw, SS, 3.5 x 10 mm, Strl
AR-8835VCL-12	VAL KreuLock™ Comp. Screw, SS, 3.5 x 12 mm
AR-8835VCL-12S	VAL KreuLock™ Screw, SS, 3.5 x 12 mm, Strl
AR-8835VCL-14	VAL KreuLock™ Comp. Screw, SS, 3.5 x 14 mm
AR-8835VCL-14S	VAL KreuLock™ Screw, SS, 3.5 x 14 mm, Strl
AR-8835VCL-16	VAL KreuLock™ Comp. Screw, SS, 3.5 x 16 mm
AR-8835VCL-16S	VAL KreuLock™ Screw, SS, 3.5 x 16 mm, Strl
AR-8835VCL-18	VAL KreuLock™ Comp. Screw, SS, 3.5 x 18 mm
AR-8835VCL-18S	VAL KreuLock™ Screw, SS, 3.5 x 18 mm, Strl
AR-8835VCL-20	VAL KreuLock™ Comp. Screw, SS, 3.5 x 20 mm
AR-8835VCL-20S	VAL KreuLock™ Screw, SS, 3.5 x 20 mm, Strl
AR-8835VCL-22	VAL KreuLock™ Comp. Screw, SS, 3.5 x 22 mm
AR-8835VCL-22S	VAL KreuLock™ Screw, SS, 3.5 x 22 mm, Strl
AR-8835VCL-24	VAL KreuLock™ Comp. Screw, SS, 3.5 x 24 mm
AR-8835VCL-24S	VAL KreuLock™ Screw, SS, 3.5 x 24 mm, Strl
AR-8835VCL-26	VAL KreuLock™ Comp. Screw, SS, 3.5 x 26 mm
AR-8835VCL-26S	VAL KreuLock™ Screw, SS, 3.5 x 26 mm, Strl
AR-8835VCL-28	VAL KreuLock™ Comp. Screw, SS, 3.5 x 28 mm
AR-8835VCL-28S	VAL KreuLock™ Screw, SS, 3.5 x 28 mm, Strl
AR-8835VCL-30	VAL KreuLock™ Comp. Screw, SS, 3.5 x 30 mm
AR-8835VCL-30S	VAL KreuLock™ Screw, SS, 3.5 x 30 mm, Strl
AR-8835VCL-32	VAL KreuLock™ Comp. Screw, SS, 3.5 x 32 mm
AR-8835VCL-32S	VAL KreuLock™ Screw, SS, 3.5 x 32 mm, Strl
AR-8835VCL-34	VAL KreuLock™ Comp. Screw, SS, 3.5 x 34 mm
AR-8835VCL-34S	VAL KreuLock™ Screw, SS, 3.5 x 34 mm, Strl
AR-8835VCL-36	VAL KreuLock™ Comp. Screw, SS, 3.5 x 36 mm
AR-8835VCL-36S	VAL KreuLock™ Screw, SS, 3.5 x 36 mm, Strl
AR-8835VCL-38	VAL KreuLock™ Comp. Screw, SS, 3.5 x 38 mm
AR-8835VCL-38S	VAL KreuLock™ Screw, SS, 3.5 x 38 mm, Strl
AR-8835VCL-40	VAL KreuLock™ Comp. Screw, SS, 3.5 x 40 mm
AR-8835VCL-40S	VAL KreuLock™ Screw, SS, 3.5 x 40 mm, Strl
AR-8835VCL-42	VAL KreuLock™ Comp. Screw, SS, 3.5 x 42 mm
AR-8835VCL-42S	VAL KreuLock™ Screw, SS, 3.5 x 42 mm, Strl

Item Number	Product Description
AR-8835VCL-44	VAL KreuLock™ Comp. Screw, SS, 3.5 x 44 mm
AR-8835VCL-44S	VAL KreuLock™ Screw, SS, 3.5 x 44 mm, Strl
AR-8835VCL-45	VAL KreuLock™ Comp. Screw, SS, 3.5 x 45 mm
AR-8835VCL-45S	VAL KreuLock™ Screw, SS, 3.5 x 45 mm, Strl
AR-8835VCL-46	VAL KreuLock™ Comp. Screw, SS, 3.5 x 46 mm
AR-8835VCL-46S	VAL KreuLock™ Screw, SS, 3.5 x 46 mm, Strl
AR-8835VCL-48	VAL KreuLock™ Comp. Screw, SS, 3.5 x 48 mm
AR-8835VCL-48S	VAL KreuLock™ Screw, SS, 3.5 x 48 mm, Strl
AR-8835VCL-50	VAL KreuLock™ Comp. Screw, SS, 3.5 x 50 mm
AR-8835VCL-50S	VAL KreuLock™ Screw, SS, 3.5 x 50 mm, Strl
AR-8835VCL-52	VAL KreuLock™ Comp. Screw, SS, 3.5 x 52 mm
AR-8835VCL-52S	VAL KreuLock™ Screw, SS, 3.5 x 52 mm, Strl
AR-8835VCL-54	VAL KreuLock™ Comp. Screw, SS, 3.5 x 54 mm
AR-8835VCL-54S	VAL KreuLock™ Screw, SS, 3.5 x 54 mm, Strl
AR-8835VCL-55	VAL KreuLock™ Comp. Screw, SS, 3.5 x 55 mm
AR-8835VCL-55S	VAL KreuLock™ Screw, SS, 3.5 x 55 mm, Strl
AR-8835VCL-56	VAL KreuLock™ Comp. Screw, SS, 3.5 x 56 mm
AR-8835VCL-56S	VAL KreuLock™ Screw, SS, 3.5 x 56 mm, Strl
AR-8835VCL-58	VAL KreuLock™ Comp. Screw, SS, 3.5 x 58 mm
AR-8835VCL-58S	VAL KreuLock™ Screw, SS, 3.5 x 58 mm, Strl
AR-8835VCL-60	VAL KreuLock™ Comp. Screw, SS, 3.5 x 60 mm
AR-8835VCL-60S	VAL KreuLock™ Screw, SS, 3.5 x 60 mm, Strl
AR-8835VCL-62	VAL KreuLock™ Comp. Screw, SS, 3.5 x 62 mm
AR-8835VCL-62S	VAL KreuLock™ Screw, SS, 3.5 x 62 mm, Strl
AR-8835VCL-64	VAL KreuLock™ Comp. Screw, SS, 3.5 x 64 mm
AR-8835VCL-64S	VAL KreuLock™ Screw, SS, 3.5 x 64 mm, Strl
AR-8835VCL-65	VAL KreuLock™ Comp. Screw, SS, 3.5 x 65 mm
AR-8835VCL-65S	VAL KreuLock™ Screw, SS, 3.5 x 65 mm, Strl
AR-8835VCL-66	VAL KreuLock™ Comp. Screw, SS, 3.5 x 66 mm
AR-8835VCL-66S	VAL KreuLock™ Screw, SS, 3.5 x 66 mm, Strl
AR-8835VCL-68	VAL KreuLock™ Comp. Screw, SS, 3.5 x 68 mm
AR-8835VCL-68S	VAL KreuLock™ Screw, SS, 3.5 x 68 mm, Strl
AR-8835VCL-70	VAL KreuLock™ Comp. Screw, SS, 3.5 x 70 mm
AR-8835VCL-70S	VAL KreuLock™ Screw, SS, 3.5 x 70 mm, Strl
AR-8835VCL-72	VAL KreuLock™ Comp. Screw, SS, 3.5 x 72 mm
AR-8835VCL-72S	VAL KreuLock™ Screw, SS, 3.5 x 72 mm, Strl
AR-8835VCL-74	VAL KreuLock™ Comp. Screw, SS, 3.5 x 74 mm
AR-8835VCL-74S	VAL KreuLock™ Screw, SS, 3.5 x 74 mm, Strl
AR-8835VCL-75	VAL KreuLock™ Comp. Screw, SS, 3.5 x 75 mm
AR-8835VCL-75S	VAL KreuLock™ Screw, SS, 3.5 x 75 mm, Strl

Item Number	Product Description
AR-8835VCL-76	VAL KreuLock™ Comp. Screw, SS, 3.5 x 76 mm
AR-8835VCL-76S	VAL KreuLock™ Screw, SS, 3.5 x 76 mm, Strl
AR-8835VCL-78	VAL KreuLock™ Comp. Screw, SS, 3.5 x 78 mm
AR-8835VCL-78S	VAL KreuLock™ Screw, SS, 3.5 x 78 mm, Strl
AR-8835VCL-80	VAL KreuLock™ Comp. Screw, SS, 3.5 x 80 mm
AR-8835VCL-80S	VAL KreuLock™ Screw, SS, 3.5 x 80 mm, Strl
AR-8835VCL-85	VAL KreuLock™ Comp. Screw, SS, 3.5 x 85 mm
AR-8835VCL-85S	VAL KreuLock™ Screw, SS, 3.5 x 85 mm, Strl
AR-8835VCL-90	VAL KreuLock™ Comp. Screw, SS, 3.5 x 90 mm
AR-8835VCL-90S	VAL KreuLock™ Screw, SS, 3.5 x 90 mm, Strl
AR-8835VCL-95	VAL KreuLock™ Comp. Screw, SS, 3.5 x 95 mm
AR-8835VCL-95S	VAL KreuLock™ Screw, SS, 3.5 x 95 mm, Strl
AR-8835VCL-100	VAL KreuLock™ Comp. Screw, SS, 3.5 x 100 mm
AR-8835VCL-100S	VAL KreuLock™ Screw, SS, 3.5 x 100 mm, Strl
AR-8835VCL-105	VAL KreuLock™ Comp. Screw, SS, 3.5 x 105 mm
AR-8835VCL-105S	VAL KreuLock™ Screw, SS, 3.5 x 105 mm, Strl
AR-8835VCL-110	VAL KreuLock™ Comp. Screw, SS, 3.5 x 110 mm
AR-8835VCL-110S	VAL KreuLock™ Screw, SS, 3.5 x 110 mm, Strl

Implant Models

VAL KreuLock™ Screw, 3.5 mm, Titanium (Ti) series, Reinforced (Ri)

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8935RVCL-10	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 10 mm, Ri
AR-8935RVCL-10S	VAL KreuLock™ Screw, Ti, 3.5 x 10 mm, Ri, Strl
AR-8935RVCL-12	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 12 mm, Ri
AR-8935RVCL-12S	VAL KreuLock™ Screw, Ti, 3.5 x 12 mm, Ri, Strl
AR-8935RVCL-14	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 14 mm, Ri
AR-8935RVCL-14S	VAL KreuLock™ Screw, Ti, 3.5 x 14 mm, Ri, Strl
AR-8935RVCL-16	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 16 mm, Ri
AR-8935RVCL-16S	VAL KreuLock™ Screw, Ti, 3.5 x 16 mm, Ri, Strl
AR-8935RVCL-18	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 18 mm, Ri
AR-8935RVCL-18S	VAL KreuLock™ Screw, Ti, 3.5 x 18 mm, Ri, Strl
AR-8935RVCL-20	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 20 mm, Ri
AR-8935RVCL-20S	VAL KreuLock™ Screw, Ti, 3.5 x 20 mm, Ri, Strl
AR-8935RVCL-22	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 22 mm, Ri
AR-8935RVCL-22S	VAL KreuLock™ Screw, Ti, 3.5 x 22 mm, Ri, Strl
AR-8935RVCL-24	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 24 mm, Ri
AR-8935RVCL-24S	VAL KreuLock™ Screw, Ti, 3.5 x 24 mm, Ri, Strl
AR-8935RVCL-26	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 26 mm, Ri
AR-8935RVCL-26S	VAL KreuLock™ Screw, Ti, 3.5 x 26 mm, Ri, Strl
AR-8935RVCL-28	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 28 mm, Ri
AR-8935RVCL-28S	VAL KreuLock™ Screw, Ti, 3.5 x 28 mm, Ri, Strl
AR-8935RVCL-30	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 30 mm, Ri
AR-8935RVCL-30S	VAL KreuLock™ Screw, Ti, 3.5 x 30 mm, Ri, Strl
AR-8935RVCL-32	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 32 mm, Ri
AR-8935RVCL-32S	VAL KreuLock™ Screw, Ti, 3.5 x 32 mm, Ri, Strl
AR-8935RVCL-34	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 34 mm, Ri
AR-8935RVCL-34S	VAL KreuLock™ Screw, Ti, 3.5 x 34 mm, Ri, Strl
AR-8935RVCL-36	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 36 mm, Ri
AR-8935RVCL-36S	VAL KreuLock™ Screw, Ti, 3.5 x 36 mm, Ri, Strl
AR-8935RVCL-38	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 38 mm, Ri
AR-8935RVCL-38S	VAL KreuLock™ Screw, Ti, 3.5 x 38 mm, Ri, Strl
AR-8935RVCL-40	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 40 mm, Ri
AR-8935RVCL-40S	VAL KreuLock™ Screw, Ti, 3.5 x 40 mm, Ri, Strl
AR-8935RVCL-42	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 42 mm, Ri
AR-8935RVCL-42S	VAL KreuLock™ Screw, Ti, 3.5 x 42 mm, Ri, Strl

Item Number	Product Description
AR-8935RVCL-44	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 44 mm, Ri
AR-8935RVCL-44S	VAL KreuLock™ Screw, Ti, 3.5 x 44 mm, RI, Strl
AR-8935RVCL-45	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 45 mm, Ri
AR-8935RVCL-45S	VAL KreuLock™ Screw, Ti, 3.5 x 45 mm, RI, Strl
AR-8935RVCL-46	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 46 mm, Ri
AR-8935RVCL-46S	VAL KreuLock™ Screw, Ti, 3.5 x 46 mm, RI, Strl
AR-8935RVCL-48	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 48 mm, Ri
AR-8935RVCL-48S	VAL KreuLock™ Screw, Ti, 3.5 x 48 mm, RI, Strl
AR-8935RVCL-50	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 50 mm, Ri
AR-8935RVCL-50S	VAL KreuLock™ Screw, Ti, 3.5 x 50 mm, RI, Strl
AR-8935RVCL-52	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 52 mm, Ri
AR-8935RVCL-52S	VAL KreuLock™ Screw, Ti, 3.5 x 52 mm, RI, Strl
AR-8935RVCL-54	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 54 mm, Ri
AR-8935RVCL-54S	VAL KreuLock™ Screw, Ti, 3.5 x 54 mm, RI, Strl
AR-8935RVCL-55	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 55 mm, Ri
AR-8935RVCL-55S	VAL KreuLock™ Screw, Ti, 3.5 x 55 mm, RI, Strl
AR-8935RVCL-56	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 56 mm, Ri
AR-8935RVCL-56S	VAL KreuLock™ Screw, Ti, 3.5 x 56 mm, RI, Strl
AR-8935RVCL-58	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 58 mm, Ri
AR-8935RVCL-58S	VAL KreuLock™ Screw, Ti, 3.5 x 58 mm, RI, Strl
AR-8935RVCL-60	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 60 mm, Ri
AR-8935RVCL-60S	VAL KreuLock™ Screw, Ti, 3.5 x 60 mm, RI, Strl
AR-8935RVCL-62	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 62 mm, Ri
AR-8935RVCL-62S	VAL KreuLock™ Screw, Ti, 3.5 x 62 mm, RI, Strl
AR-8935RVCL-64	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 64 mm, Ri
AR-8935RVCL-64S	VAL KreuLock™ Screw, Ti, 3.5 x 64 mm, RI, Strl
AR-8935RVCL-65	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 65 mm, Ri
AR-8935RVCL-65S	VAL KreuLock™ Screw, Ti, 3.5 x 65 mm, RI, Strl
AR-8935RVCL-66	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 66 mm, Ri
AR-8935RVCL-66S	VAL KreuLock™ Screw, Ti, 3.5 x 66 mm, RI, Strl
AR-8935RVCL-68	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 68 mm, Ri
AR-8935RVCL-68S	VAL KreuLock™ Screw, Ti, 3.5 x 68 mm, RI, Strl
AR-8935RVCL-70	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 70 mm, Ri
AR-8935RVCL-70S	VAL KreuLock™ Screw, Ti, 3.5 x 70 mm, RI, Strl
AR-8935RVCL-72	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 72 mm, Ri
AR-8935RVCL-72S	VAL KreuLock™ Screw, Ti, 3.5 x 72 mm, RI, Strl
AR-8935RVCL-74	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 74 mm, Ri
AR-8935RVCL-74S	VAL KreuLock™ Screw, Ti, 3.5 x 74 mm, RI, Strl
AR-8935RVCL-75	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 75 mm, Ri
AR-8935RVCL-75S	VAL KreuLock™ Screw, Ti, 3.5 x 75 mm, RI, Strl

Item Number	Product Description
AR-8935RVCL-76	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 76 mm, Ri
AR-8935RVCL-76S	VAL KreuLock™ Screw, Ti, 3.5 x 76 mm, Ri, Strl
AR-8935RVCL-78	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 78 mm, Ri
AR-8935RVCL-78S	VAL KreuLock™ Screw, Ti, 3.5 x 78 mm, Ri, Strl
AR-8935RVCL-80	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 80 mm, Ri
AR-8935RVCL-80S	VAL KreuLock™ Screw, Ti, 3.5 x 80 mm, Ri, Strl
AR-8935RVCL-85	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 85 mm, Ri
AR-8935RVCL-85S	VAL KreuLock™ Screw, Ti, 3.5 x 85 mm, Ri, Strl
AR-8935RVCL-90	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 90 mm, Ri
AR-8935RVCL-90S	VAL KreuLock™ Screw, Ti, 3.5 x 90 mm, Ri, Strl
AR-8935RVCL-95	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 95 mm, Ri
AR-8935RVCL-95S	VAL KreuLock™ Screw, Ti, 3.5 x 95 mm, Ri, Strl
AR-8935RVCL-100	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 100 mm, Ri
AR-8935RVCL-100S	VAL KreuLock™ Screw, Ti, 3.5 x 100 mm, Ri, Strl
AR-8935RVCL-105	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 105 mm, Ri
AR-8935RVCL-105S	VAL KreuLock™ Screw, Ti, 3.5 x 105 mm, Ri, Strl
AR-8935RVCL-110	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 110 mm, Ri
AR-8935RVCL-110S	VAL KreuLock™ Screw, Ti, 3.5 x 110 mm, Ri, Strl

Implant Models

VAL KreuLock™ Screw, 3.5 mm, Titanium (Ti) series

Consult your VAL KreuLock Screw Patient Implant Card (PIC) for information on the device type/model of the implant used in your procedure.

Item Number	Product Description
AR-8935VCL-10	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 10 mm
AR-8935VCL-10S	VAL KreuLock™ Screw, Ti, 3.5 x 10 mm, Strl
AR-8935VCL-12	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 12 mm
AR-8935VCL-12S	VAL KreuLock™ Screw, Ti, 3.5 x 12 mm, Strl
AR-8935VCL-14	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 14 mm
AR-8935VCL-14S	VAL KreuLock™ Screw, Ti, 3.5 x 14 mm, Strl
AR-8935VCL-16	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 16 mm
AR-8935VCL-16S	VAL KreuLock™ Screw, Ti, 3.5 x 16 mm, Strl
AR-8935VCL-18	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 18 mm
AR-8935VCL-18S	VAL KreuLock™ Screw, Ti, 3.5 x 18 mm, Strl
AR-8935VCL-20	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 20 mm
AR-8935VCL-20S	VAL KreuLock™ Screw, Ti, 3.5 x 20 mm, Strl
AR-8935VCL-22	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 22 mm
AR-8935VCL-22S	VAL KreuLock™ Screw, Ti, 3.5 x 22 mm, Strl
AR-8935VCL-24	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 24 mm
AR-8935VCL-24S	VAL KreuLock™ Screw, Ti, 3.5 x 24 mm, Strl
AR-8935VCL-26	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 26 mm
AR-8935VCL-26S	VAL KreuLock™ Screw, Ti, 3.5 x 26 mm, Strl
AR-8935VCL-28	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 28 mm
AR-8935VCL-28S	VAL KreuLock™ Screw, Ti, 3.5 x 28 mm, Strl
AR-8935VCL-30	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 30 mm
AR-8935VCL-30S	VAL KreuLock™ Screw, Ti, 3.5 x 30 mm, Strl
AR-8935VCL-32	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 32 mm
AR-8935VCL-32S	VAL KreuLock™ Screw, Ti, 3.5 x 32 mm, Strl
AR-8935VCL-34	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 34 mm
AR-8935VCL-34S	VAL KreuLock™ Screw, Ti, 3.5 x 34 mm, Strl
AR-8935VCL-36	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 36 mm
AR-8935VCL-36S	VAL KreuLock™ Screw, Ti, 3.5 x 36 mm, Strl
AR-8935VCL-38	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 38 mm
AR-8935VCL-38S	VAL KreuLock™ Screw, Ti, 3.5 x 38 mm, Strl
AR-8935VCL-40	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 40 mm
AR-8935VCL-40S	VAL KreuLock™ Screw, Ti, 3.5 x 40 mm, Strl
AR-8935VCL-42	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 42 mm
AR-8935VCL-42S	VAL KreuLock™ Screw, Ti, 3.5 x 42 mm, Strl

Item Number	Product Description
AR-8935VCL-44	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 44 mm
AR-8935VCL-44S	VAL KreuLock™ Screw, Ti, 3.5 x 44 mm, Strl
AR-8935VCL-45	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 45 mm
AR-8935VCL-45S	VAL KreuLock™ Screw, Ti, 3.5 x 45 mm, Strl
AR-8935VCL-46	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 46 mm
AR-8935VCL-46S	VAL KreuLock™ Screw, Ti, 3.5 x 46 mm, Strl
AR-8935VCL-48	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 48 mm
AR-8935VCL-48S	VAL KreuLock™ Screw, Ti, 3.5 x 48 mm, Strl
AR-8935VCL-50	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 50 mm
AR-8935VCL-50S	VAL KreuLock™ Screw, Ti, 3.5 x 50 mm, Strl
AR-8935VCL-52	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 52 mm
AR-8935VCL-52S	VAL KreuLock™ Screw, Ti, 3.5 x 52 mm, Strl
AR-8935VCL-54	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 54 mm
AR-8935VCL-54S	VAL KreuLock™ Screw, Ti, 3.5 x 54 mm, Strl
AR-8935VCL-55	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 55 mm
AR-8935VCL-55S	VAL KreuLock™ Screw, Ti, 3.5 x 55 mm, Strl
AR-8935VCL-56	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 56 mm
AR-8935VCL-56S	VAL KreuLock™ Screw, Ti, 3.5 x 56 mm, Strl
AR-8935VCL-58	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 58 mm
AR-8935VCL-58S	VAL KreuLock™ Screw, Ti, 3.5 x 58 mm, Strl
AR-8935VCL-60	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 60 mm
AR-8935VCL-60S	VAL KreuLock™ Screw, Ti, 3.5 x 60 mm, Strl
AR-8935VCL-62	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 62 mm
AR-8935VCL-62S	VAL KreuLock™ Screw, Ti, 3.5 x 62 mm, Strl
AR-8935VCL-64	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 64 mm
AR-8935VCL-64S	VAL KreuLock™ Screw, Ti, 3.5 x 64 mm, Strl
AR-8935VCL-65	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 65 mm
AR-8935VCL-65S	VAL KreuLock™ Screw, Ti, 3.5 x 65 mm, Strl
AR-8935VCL-66	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 66 mm
AR-8935VCL-66S	VAL KreuLock™ Screw, Ti, 3.5 x 66 mm, Strl
AR-8935VCL-68	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 68 mm
AR-8935VCL-68S	VAL KreuLock™ Screw, Ti, 3.5 x 68 mm, Strl
AR-8935VCL-70	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 70 mm
AR-8935VCL-70S	VAL KreuLock™ Screw, Ti, 3.5 x 70 mm, Strl
AR-8935VCL-72	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 72 mm
AR-8935VCL-72S	VAL KreuLock™ Screw, Ti, 3.5 x 72 mm, Strl
AR-8935VCL-74	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 74 mm
AR-8935VCL-74S	VAL KreuLock™ Screw, Ti, 3.5 x 74 mm, Strl
AR-8935VCL-75	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 75 mm
AR-8935VCL-75S	VAL KreuLock™ Screw, Ti, 3.5 x 75 mm, Strl

Item Number	Product Description
AR-8935VCL-76	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 76 mm
AR-8935VCL-76S	VAL KreuLock™ Screw, Ti, 3.5 x 76 mm, Strl
AR-8935VCL-78	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 78 mm
AR-8935VCL-78S	VAL KreuLock™ Screw, Ti, 3.5 x 78 mm, Strl
AR-8935VCL-80	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 80 mm
AR-8935VCL-80S	VAL KreuLock™ Screw, Ti, 3.5 x 80 mm, Strl
AR-8935VCL-85	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 85 mm
AR-8935VCL-85S	VAL KreuLock™ Screw, Ti, 3.5 x 85 mm, Strl
AR-8935VCL-90	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 90 mm
AR-8935VCL-90S	VAL KreuLock™ Screw, Ti, 3.5 x 90 mm, Strl
AR-8935VCL-95	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 95 mm
AR-8935VCL-95S	VAL KreuLock™ Screw, Ti, 3.5 x 95 mm, Strl
AR-8935VCL-100	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 100 mm
AR-8935VCL-100S	VAL KreuLock™ Screw, Ti, 3.5 x 100 mm, Strl
AR-8935VCL-105	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 105 mm
AR-8935VCL-105S	VAL KreuLock™ Screw, Ti, 3.5 x 105 mm, Strl
AR-8935VCL-110	VAL KreuLock™ Comp. Screw, Ti, 3.5 x 110 mm
AR-8935VCL-110S	VAL KreuLock™ Screw, Ti, 3.5 x 110 mm, Strl

Contact Information

Any serious incident that occurs in relation to the device should be reported to the manufacturer and to the health authority where the incident occurred.

Region	Contact
 Arthrex, Inc.	1370 Creekside Blvd. Naples, FL 34108 USA ☎ +1 800 934-4404 www.arthrex.com
Arthrex GmbH	Erwin-Hielscher-Strasse 9 81249 München Germany ☎ +49 89 90 90 05-0 www.arthrex.de

USA – U. S. Food & Drug Administration website: <https://www.fda.gov/medical-devices/resources-you-medical-devices/consumers-medical-devices>

Symbols glossary can be found at www.arthrex.com/symbolsglossary.



The information contained in this patient information leaflet is not medical advice and is not meant to be a substitute for the advice provided by a surgeon or other qualified medical professional on the use of these products. You should talk with your physician or healthcare provider for more information about your health condition, and whether Arthrex products might be appropriate for you. The surgeon who performs any surgical procedure is responsible for determining and using the appropriate techniques for surgical procedures on each individual patient. Arthrex recommends that surgeons be trained on the use of any particular product before using it in surgery. A surgeon must always rely on their professional medical judgment when deciding whether to use a particular product when treating a particular patient. A surgeon must always refer to the package insert, product label, and/or directions for use before using any Arthrex product. Products may not be available in all markets because product availability is subject to regulatory approvals and medical practices in individual markets. Please contact Arthrex if you have questions about the availability of products in your area.

arthrex.com

© 2026-08 Arthrex, Inc. All rights reserved. PIL-0019_USr2-en