

Fluid Diamond^{®*} Floor Aspirator

Contain Fluids in the OR



The Fluid Diamond disposable floor aspirator is designed to quickly and quietly remove arthroscopic fluids from the operating room to support a safe working environment. With a unique shape, the Fluid Diamond floor aspirator can slide further underneath a stationary object to remove hard-to-reach fluids and reach into corners where round designs cannot. Multiple ports catch fluids, which are quickly channeled to the bottom of the Fluid Diamond aspirator and evacuated. The Fluid Diamond aspirator comes with suction tubing and can be attached to a vacuum suction, a portable suction device, or directly to The Suctioner[®] or The EcoSuctioner[®] surgical suction mats.

- Fast, effective fluid collection
- Powerful suction action keeps floors dry and reduces the potential for slipping or falling in the OR
- Attaches to any surgical suction source
- Eliminates the need for slippery, costly towels and blankets on the OR floor to absorb fluids
- Includes suction tubing for easy attachment
- Designed with quieting pad and anti-floor-sticking technology

Product Description	Item Number
Fluid Diamond Floor Aspirator, qty. 30	USA-FDX

*The Fluid Diamond, The Suctioner, and The EcoSuctioner are trademarks of Flagship Surgical, LLC.



arthrex.com

© 2021 Arthrex, Inc. All rights reserved. LS1-00039-en-US_C

Arthrex[®] 