



**COOLCUT RADIO FREQUENCY ABLATION**

**CoolCut™**  
**RADIO FREQUENCY ABLATION**


# Arthrex CoolCut Radio Frequency Ablation

The new line of CoolCut radio frequency devices features an innovative ergonomic handle with finger touch controls. The 90° versions have a unique shape that allows for maximum extension from the tip of the electrode to the shaft of the device, while still fitting through a 5 mm cannula. The new CoolCut shaft coating facilitates easy insertion through both working cannulas and subcutaneous tissue.

The pioneering design of the CoolCut electrode tips makes it possible to operate the CoolCut series at extremely low power settings, thereby increasing patient safety without sacrificing performance. The reduction of the power requirements also allows for the absence of a ceramic ring surrounding the electrode face, reducing the size of the tip which aids in reaching target tissue in anatomically tight spaces.

## **Arthrex CoolCut**

The Next Generation of Radio Frequency Ablation Devices is Here Now.



## CoolCut Ablation Wands

### Aspirating Ablation Probes



**CoolCut 90**  
AR-9803A-90



**CoolCut 50**  
AR-9803A-50



**CoolCut 30**  
AR-9803A-30

### Meniscal Probes



**CoolCut Meniscectomy 45**  
AR-9802M-45



**CoolCut Meniscectomy 90**  
AR-9802M-90

### Hook Tip Probes



**CoolCut Hook**  
AR-9801

### Non-Aspirating Ablation Probes



**CoolCut 90NA**  
AR-9803-90



**CoolCut 50NA**  
AR-9803-50



**CoolCut 30NA**  
AR-9803-30

### Ball Probes



**CoolCut Tip**  
AR-9808



**CoolCut Tip 45**  
AR-9808-45

## Accessory Devices

AR-9610

AR-9610SGP

Arthrex Pencil Kit

Split Ground Pad

## OPES Ablation Wands



**ToothBrush**  
AR-9705A-90



**ToothBrush NA**  
AR-9705-90



**3mm Low Profile**  
AR-9703A-90

## Small Joint CoolCut Ablation Wands



**CoolCut SJ Ball**  
AR-9808SJ



**CoolCut SJ Ball 45**  
AR-9808SJ-45



**CoolCut SJ 45**  
AR-9809SJ-45



**CoolCut SJ 90**  
AR-9809SJ-90



©2013, Arthrex Inc. All rights reserved.  
LB1-0670-EN\_B

***<http://raisingthebar.Arthrex.com>***